




OR4C13 Antibody

Product Code	CSB-PA003548
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGP0
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 4C13.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR4C13; Olfactory receptor 4C13; Olfactory receptor OR11-260
Immunogen Species	Homo sapiens (Human)
Target Names	OR4C13

Image



Western Blot analysis of MCF7 cells using Olfactory receptor 4C13 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.