




OR4C6 Antibody

Product Code	CSB-PA003550
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NH72
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 4C6.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR4C6; Olfactory receptor 4C6; Olfactory receptor OR11-138
Immunogen Species	Homo sapiens (Human)
Target Names	OR4C6

Image



Western Blot analysis of HepG2 cells using Olfactory receptor 4C6 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.