




OR4Q3 Antibody

Product Code	CSB-PA003555
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NH05
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 4Q3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR4Q3; C14orf13; OR4Q4; Olfactory receptor 4Q3; Olfactory receptor 4Q4; Olfactory receptor OR14-3
Immunogen Species	Homo sapiens (Human)
Target Names	OR4Q3

Image



Western Blot analysis of HepG2 cells using Olfactory receptor 4Q3 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.