




OR4Q3 Antibody

Product Code	CSB-PA003555
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NH05
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 4Q3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR4Q3; C14orf13; OR4Q4; Olfactory receptor 4Q3; Olfactory receptor 4Q4; Olfactory receptor OR14-3
Immunogen Species	Homo sapiens (Human)
Target Names	OR4Q3

Image



Western Blot analysis of HepG2 cells using Olfactory receptor 4Q3 Polyclonal Antibody