




OR52B2 Antibody

Product Code	CSB-PA003562
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RD2
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 52B2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR52B2; Olfactory receptor 52B2; Olfactory receptor OR11-70
Immunogen Species	Homo sapiens (Human)
Target Names	OR52B2

Image



Western Blot analysis of HepG2 cells using Olfactory receptor 52B2 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.