




OR52W1 Antibody

Product Code	CSB-PA003570
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IF63
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 52W1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR52W1; OR52W1P; Olfactory receptor 52W1; Olfactory receptor OR11-71
Immunogen Species	Homo sapiens (Human)
Target Names	OR52W1

Image



Western Blot analysis of HeLa cells using Olfactory receptor 52W1 Polyclonal Antibody