




OR56B4 Antibody

Product Code	CSB-PA003571
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NH76
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 56B4.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR56B4; Olfactory receptor 56B4; Olfactory receptor OR11-67
Immunogen Species	Homo sapiens (Human)
Target Names	OR56B4

Image



Western Blot analysis of COLO205 cells using Olfactory receptor 56B4 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.