




OR5AK3P Antibody

Product Code	CSB-PA003572
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NH89
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 5AK3.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR5AK3P; OR5AK3; Putative olfactory receptor 5AK3
Immunogen Species	Homo sapiens (Human)
Target Names	OR5AK3P

Image



Western Blot analysis of COS7 cells using Olfactory receptor 5AK3 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.