




OR5P3 Antibody

Product Code	CSB-PA003577
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WZ94
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 5P3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR5P3; Olfactory receptor 5P3; Olfactory receptor OR11-94; Olfactory receptor-like protein JCG1
Immunogen Species	Homo sapiens (Human)
Target Names	OR5P3

Image



Western Blot analysis of COLO205 cells using Olfactory receptor 5P3 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.