




OR6C68 Antibody

Product Code	CSB-PA003579
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NDL8
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 6C68.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR6C68; Olfactory receptor 6C68
Immunogen Species	Homo sapiens (Human)
Target Names	OR6C68

Image



Western Blot analysis of HuvEc cells using Olfactory receptor 6C68 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.