




OR7A10 Antibody

Product Code	CSB-PA003582
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76100
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 7A10.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR7A10; Olfactory receptor 7A10; OST027; Olfactory receptor OR19-18
Immunogen Species	Homo sapiens (Human)
Target Names	OR7A10

Image



Western Blot analysis of COS7 cells using Olfactory receptor 7A10 Polyclonal Antibody