




OR7C1 Antibody

Product Code	CSB-PA003583
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76099
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 7C1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR7C1; OR7C4; Olfactory receptor 7C1; Olfactory receptor 7C4; Olfactory receptor OR19-16; Olfactory receptor TPCR86
Immunogen Species	Homo sapiens (Human)
Target Names	OR7C1

Image



Western Blot analysis of COLO205 cells using Olfactory receptor 7C1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.