




OR7C2 Antibody

Product Code	CSB-PA003584
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60412
Immunogen	Synthesized peptide derived from the Internal region of Human Olfactory receptor 7C2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR7C2; OR7C3; Olfactory receptor 7C2; Olfactory receptor 19-18; OR19-18; Olfactory receptor 7C3; Olfactory receptor OR19-22
Immunogen Species	Homo sapiens (Human)
Target Names	OR7C2

Image



Western Blot analysis of HuvEc cells using Olfactory receptor 7C2 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.