




OR8J3 Antibody

Product Code	CSB-PA003587
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NGG0
Immunogen	Synthesized peptide derived from the C-terminal region of Human Olfactory receptor 8J3.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	OR8J3; Olfactory receptor 8J3; Olfactory receptor OR11-173
Immunogen Species	Homo sapiens (Human)
Target Names	OR8J3

Image



Western Blot analysis of HuvEc cells using Olfactory receptor 8J3 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.