




TNFRSF11B Antibody

Product Code	CSB-PA003597
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00300
Immunogen	Synthesized peptide derived from the N-terminal region of Human OPG.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TNFRSF11B; OCIF; OPG; Tumor necrosis factor receptor superfamily member 11B; Osteoclastogenesis inhibitory factor; Osteoprotegerin
Immunogen Species	Homo sapiens (Human)
Target Names	TNFRSF11B

Image



Western Blot analysis of 293 cells using OPG Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.