





USO1 Antibody

Product Code	CSB-PA003610
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60763
Immunogen	Synthesized peptide derived from the C-terminal region of Human p115.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	USO1; VDP; General vesicular transport factor p115; Protein USO1 homolog; Transcytosis-associated protein; TAP; Vesicle-docking protein
Immunogen Species	Homo sapiens (Human)
Target Names	USO1

Image



Western Blot analysis of K562 cells using p115 Polyclonal Antibody



Western Blot analysis of K562 293 cells using p115 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.