




BCAR1 Antibody

Product Code	CSB-PA003613
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56945
Immunogen	Synthesized peptide derived from Human p130 Cas around the non-phosphorylation site of Y165.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BCAR1; CAS; CASS1; CRKAS; Breast cancer anti-estrogen resistance protein 1; CRK-associated substrate; Cas scaffolding protein family member 1; p130cas
Immunogen Species	Homo sapiens (Human)
Target Names	BCAR1

Image



Western Blot analysis of HepG2 cells using p130 Cas Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.