




CDKN2C Antibody

Product Code	CSB-PA003620
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42773
Immunogen	Synthesized peptide derived from the C-terminal region of Human p18.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CDKN2C; CDKN6; Cyclin-dependent kinase 4 inhibitor C; Cyclin-dependent kinase 6 inhibitor; p18-INK4c; p18-INK6
Immunogen Species	Homo sapiens (Human)
Target Names	CDKN2C

Image



Western Blot analysis of COS7 cells using p18 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.