





CDKN1A Antibody

Product Code	CSB-PA003623
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38936
Immunogen	Synthesized peptide derived from Human p21 around the non-phosphorylation site of T145.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CDKN1A; CAP20; CDKN1; CIP1; MDA6; PIC1; SDI1; WAF1; Cyclin-dependent kinase inhibitor 1; CDK-interacting protein 1; Melanoma differentiation-associated protein 6; MDA-6; p21
Immunogen Species	Homo sapiens (Human)
Target Names	CDKN1A

Image



Western Blot analysis of COLO205 cells using p21 Polyclonal Antibody



Western Blot analysis of HepG2 HELA cells using p21 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.