





CDKN1B Antibody

Product Code	CSB-PA003627
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46527
Immunogen	Synthesized peptide derived from Human p27 around the non-phosphorylation site of T198.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CDKN1B; KIP1; Cyclin-dependent kinase inhibitor 1B; Cyclin-dependent kinase inhibitor p27; p27Kip1
Immunogen Species	Homo sapiens (Human)
Target Names	CDKN1B

Image



Western Blot analysis of COLO205 cells using p27 Polyclonal Antibody



Western Blot analysis of KB cells using p27 Polyclonal Antibody