





CDK5R1 Antibody

Product Code	CSB-PA003660
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15078
Immunogen	Synthesized peptide derived from the N-terminal region of Human p35.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CDK5R1; CDK5R; NCK5A; Cyclin-dependent kinase 5 activator 1; CDK5 activator 1; Cyclin-dependent kinase 5 regulatory subunit 1; TPKII regulatory subunit
Immunogen Species	Homo sapiens (Human)
Target Names	CDK5R1

Image



Western Blot analysis of HELA cells using p35 Polyclonal Antibody



Western Blot analysis of K562 cells using p35 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.