




CDK5R2 Antibody

Product Code	CSB-PA003663
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13319
Immunogen	Synthesized peptide derived from the Internal region of Human p39.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CDK5R2; NCK5A1; Cyclin-dependent kinase 5 activator 2; CDK5 activator 2; Cyclin-dependent kinase 5 regulatory subunit 2; p39; p39l
Immunogen Species	Homo sapiens (Human)
Target Names	CDK5R2

Image



Western Blot analysis of HuvEc cells using p39 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.