




NCF1 Antibody

Product Code	CSB-PA003664
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14598
Immunogen	Synthesized peptide derived from Human p47-phox around the non-phosphorylation site of S304.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NCF1; NOXO2; SH3PXD1A; Neutrophil cytosol factor 1; NCF-1; 47 kDa autosomal chronic granulomatous disease protein; 47 kDa neutrophil oxidase factor; NCF-47K; Neutrophil NADPH oxidase factor 1; Nox organizer 2; Nox-organizing protein 2; SH3
Immunogen Species	Homo sapiens (Human)
Target Names	NCF1

Image



Western Blot analysis of COLO205 cells using p47-phox Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.