





# NCF1 Antibody

<b>Product Code</b>	CSB-PA003665
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P14598
<b>Immunogen</b>	Synthesized peptide derived from Human p47-phox around the non-phosphorylation site of S370.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	NCF1; NOXO2; SH3PXD1A; Neutrophil cytosol factor 1; NCF-1; 47 kDa autosomal chronic granulomatous disease protein; 47 kDa neutrophil oxidase factor; NCF-47K; Neutrophil NADPH oxidase factor 1; Nox organizer 2; Nox-organizing protein 2; SH3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	NCF1

## Image



Western Blot analysis of COLO205 cells using p47-phox Polyclonal Antibody



Western Blot analysis of HuvEc cells using p47-phox Polyclonal Antibody

### Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.