





TP53 Antibody

Product Code	CSB-PA003677
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04637
Immunogen	Synthesized peptide derived from Human p53 around the non-phosphorylation site of S378.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TP53; P53; Cellular tumor antigen p53; Antigen NY-CO-13; Phosphoprotein p53; Tumor suppressor p53
Immunogen Species	Homo sapiens (Human)
Target Names	TP53

Image



Western Blot analysis of SH-SY5Y cells using p53 Polyclonal Antibody



Western Blot analysis of HELA 823 293-UV
22RV1 cells using p53 Polyclonal Antibody