





RNF144B Antibody

Product Code	CSB-PA003683
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z419
Immunogen	Synthesized peptide derived from the Internal region of Human p53RFP.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	RNF144B; IBRDC2; P53RFP; E3 ubiquitin-protein ligase RNF144B; IBR domain-containing protein 2; RING finger protein 144B; p53-inducible RING finger protein
Immunogen Species	Homo sapiens (Human)
Target Names	RNF144B

Image



Western Blot analysis of HeLa cells using p53RFP Polyclonal Antibody



Western Blot analysis of RAT-MUSCLE HELA 293T AD293 cells using p53RFP Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.