




DDX5 Antibody

Product Code	CSB-PA003685
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17844
Immunogen	Synthesized peptide derived from Human p68 RNA Helicase around the non-phosphorylation site of Y593.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DDX5; G17P1; HELR; HLR1; Probable ATP-dependent RNA helicase DDX5; DEAD box protein 5; RNA helicase p68
Immunogen Species	Homo sapiens (Human)
Target Names	DDX5

Image



Western Blot analysis of HepG2 cells using p68 RNA Helicase Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.