





PAK7 Antibody

Product Code	CSB-PA003702
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P286
Immunogen	Synthesized peptide derived from the C-terminal region of Human PAK5.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PAK7; KIAA1264; PAK5; Serine/threonine-protein kinase PAK 7; p21-activated kinase 5; PAK-5; p21-activated kinase 7; PAK-7
Immunogen Species	Homo sapiens (Human)
Target Names	PAK5

Image



Western Blot analysis of 3T3 cells using PAK5 Polyclonal Antibody



Western Blot analysis of various cells using
PAK5 Polyclonal Antibody