




PAX5 Antibody

Product Code	CSB-PA003710
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02548
Immunogen	Synthesized peptide derived from the N-terminal region of Human Pax-5.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PAX5; Paired box protein Pax-5; B-cell-specific transcription factor; BSAP
Immunogen Species	Homo sapiens (Human)
Target Names	PAX5

Image



Western Blot analysis of 293 cells using Pax-5 Polyclonal Antibody