





PCNA Antibody

Product Code	CSB-PA003719
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12004
Immunogen	Synthesized peptide derived from the Internal region of Human PCNA.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PCNA; Proliferating cell nuclear antigen; PCNA; Cyclin
Immunogen Species	Homo sapiens (Human)
Target Names	PCNA

Image



Western Blot analysis of COS7 cells using PCNA Polyclonal Antibody



Western Blot analysis of HepG2 HuvEc cells using PCNA Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.