




PDCD4 Antibody

Product Code	CSB-PA003720
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53EL6
Immunogen	Synthesized peptide derived from Human Pdc4 around the non-phosphorylation site of S457.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PDCD4; H731; Programmed cell death protein 4; Neoplastic transformation inhibitor protein; Nuclear antigen H731-like; Protein 197/15a
Immunogen Species	Homo sapiens (Human)
Target Names	PDCD4

Image



Western Blot analysis of HeLa cells using Pdc4 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.