





PDGFRA Antibody

Product Code	CSB-PA003724
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16234
Immunogen	Synthesized peptide derived from the C-terminal region of Human PDGFR- α .
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PDGFRA; PDGFR2; RHEPDGFRA; Platelet-derived growth factor receptor alpha; PDGF-R-alpha; PDGFR-alpha; Alpha platelet-derived growth factor receptor; Alpha-type platelet-derived growth factor receptor; CD140 antigen-like family member A; CD14
Immunogen Species	Homo sapiens (Human)
Target Names	PDGFRA

Image



Western Blot analysis of HepG2 cells using PDGFR- α Polyclonal Antibody



Western Blot analysis of HepG2 cells using PDGFR- α Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.