




PDGFRB Antibody

Product Code	CSB-PA003726
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09619
Immunogen	Synthesized peptide derived from Human PDGFR- β around the non-phosphorylation site of Y751.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PDGFRB; PDGFR; PDGFR1; Platelet-derived growth factor receptor beta; PDGF-R-beta; PDGFR-beta; Beta platelet-derived growth factor receptor; Beta-type platelet-derived growth factor receptor; CD140 antigen-like family member B; Platelet-deri
Immunogen Species	Homo sapiens (Human)
Target Names	PDGFRB

Image



Western Blot analysis of HELA KB cells using PDGFR- β Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.