





PDPK1 Antibody

Product Code	CSB-PA003730
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15530
Immunogen	Synthesized peptide derived from Human PDK1 around the non-phosphorylation site of S241.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PDPK1; PDK1; 3-phosphoinositide-dependent protein kinase 1; hPDK1
Immunogen Species	Homo sapiens (Human)
Target Names	PDPK1

Image



Western Blot analysis of COLO205 cells using PDK1 Polyclonal Antibody



Western Blot analysis of HELA COLO 293T cells using PDK1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.