





PFKL Antibody

Product Code	CSB-PA003748
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17858
Immunogen	Synthesized peptide derived from the C-terminal region of Human PFK-B.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PFKL; 6-phosphofructokinase; liver type; Phosphofructo-1-kinase isozyme B; PFK-B; Phosphofructokinase 1; Phosphohexokinase
Immunogen Species	Homo sapiens (Human)
Target Names	PFKL

Image



Western Blot analysis of NIH-3T3 cells using PFK-B Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using PFK-B Polyclonal Antibody