




PGRMC2 Antibody

Product Code	CSB-PA003753
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15173
Immunogen	Synthesized peptide derived from the C-terminal region of Human PGRMC2.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PGRMC2; DG6; PMBP; Membrane-associated progesterone receptor component 2; Progesterone membrane-binding protein; Steroid receptor protein DG6
Immunogen Species	Homo sapiens (Human)
Target Names	PGRMC2

Image



Western Blot analysis of HeLa cells using PGRMC2 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.