





PIK3CA Antibody

Product Code	CSB-PA003761
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42336
Immunogen	Synthesized peptide derived from the Internal region of Human PI 3-kinase p110 α .
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PIK3CA; Phosphatidylinositol 4; 5-bisphosphate 3-kinase catalytic subunit alpha isoform; PI3-kinase subunit alpha; PI3K-alpha; PI3Kalpha; PtdIns-3-kinase subunit alpha; Phosphatidylinositol 4,5-bisphosphate 3-kinase 110 kDa catalytic subunit
Immunogen Species	Homo sapiens (Human)
Target Names	PIK3CA

Image



Western Blot analysis of HELA cells using PI 3-kinase p110 α Polyclonal Antibody



Western Blot analysis of HELA cells using PI 3-kinase p110 α Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.