





PIK3R2 Antibody

Product Code	CSB-PA003770
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00459
Immunogen	Synthesized peptide derived from Human PI 3-kinase p85 β around the non-phosphorylation site of Y464.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PIK3R2; Phosphatidylinositol 3-kinase regulatory subunit beta; PI3-kinase regulatory subunit beta; PI3K regulatory subunit beta; PtdIns-3-kinase regulatory subunit beta; Phosphatidylinositol 3-kinase 85 kDa regulatory subunit beta; PI3-kinase
Immunogen Species	Homo sapiens (Human)
Target Names	PIK3R2

Image



Western Blot analysis of 453 cells using PI 3-kinase p85 β Polyclonal Antibody



Western Blot analysis of 293T AD293 HELA HepG2 cells using PI 3-kinase p85 β Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.