





# PIK3R2 Antibody

<b>Product Code</b>	CSB-PA003770
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O00459
<b>Immunogen</b>	Synthesized peptide derived from Human PI 3-kinase p85 $\beta$ around the non-phosphorylation site of Y464.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	PIK3R2; Phosphatidylinositol 3-kinase regulatory subunit beta; PI3-kinase regulatory subunit beta; PI3K regulatory subunit beta; PtdIns-3-kinase regulatory subunit beta; Phosphatidylinositol 3-kinase 85 kDa regulatory subunit beta; PI3-kina
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	PIK3R2

## Image



Western Blot analysis of 453 cells using PI 3-kinase p85 $\beta$  Polyclonal Antibody



Western Blot analysis of 293T AD293 HELA  
HepG2 cells using PI 3-kinase p85 $\beta$  Polyclonal  
Antibody