





PIAS2 Antibody

Product Code	CSB-PA003781
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75928
Immunogen	Synthesized peptide derived from the N-terminal region of Human PIASx.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PIAS2; PIASX; E3 SUMO-protein ligase PIAS2; Androgen receptor-interacting protein 3; ARIP3; DAB2-interacting protein; DIP; Msx-interacting zinc finger protein; Miz1; PIAS-NY protein; Protein inhibitor of activated STAT x; Protein inhibitor
Immunogen Species	Homo sapiens (Human)
Target Names	PIAS2

Image



Western Blot analysis of COS7 cells using PIASx Polyclonal Antibody



Western Blot analysis of NIH-3T3 Jurkat K562 cells using PIASx Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.