




CLNS1A Antibody

Product Code	CSB-PA003782
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54105
Immunogen	Synthesized peptide derived from the C-terminal region of Human pICln.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CLNS1A; CLCI; ICLN; Methylosome subunit pICln; Chloride channel; nucleotide sensitive 1A; Chloride conductance regulatory protein ICln; I(Cln); Chloride ion current inducer protein; CICI; Reticulocyte pICln
Immunogen Species	Homo sapiens (Human)
Target Names	CLNS1A

Image



Western Blot analysis of Jurkat cells using pICln Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.