





PIP5K1C Antibody

Product Code	CSB-PA003788
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60331
Immunogen	Synthesized peptide derived from the Internal region of Human PIPK I γ .
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PIP5K1C; KIAA0589; Phosphatidylinositol 4-phosphate 5-kinase type-1 gamma; PIP5K1-gamma; PtdIns(4)P-5-kinase 1 gamma; Phosphatidylinositol 4-phosphate 5-kinase type I gamma; PIP5K1gamma
Immunogen Species	Homo sapiens (Human)
Target Names	PIP5K1C

Image



Western Blot analysis of A549 cells using PIPK I γ Polyclonal Antibody



Western Blot analysis of various cells using PIPK I γ Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.