





PRKAR2B Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA003794 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P31323 |
| Immunogen | Synthesized peptide derived from Human PKA II β reg around the non-phosphorylation site of S113. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | PRKAR2B; cAMP-dependent protein kinase type II-beta regulatory subunit |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | PRKAR2B |

Image



Western Blot analysis of 3T3 cells using PKA II β reg Polyclonal Antibody



Western Blot analysis of NIH-3T3 K562 HepG2 cells using PKA II β reg Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.