





PRKACA/PRKACB Antibody

Product Code	CSB-PA003798
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17612/P22694
Immunogen	Synthesized peptide derived from the N-terminal region of Human PKA α / β cat.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PRKACA; PKACA; cAMP-dependent protein kinase catalytic subunit alpha; PKA C-alpha; PRKACB; cAMP-dependent protein kinase catalytic subunit beta; PKA C-beta
Immunogen Species	Homo sapiens (Human)
Target Names	PRKACA/PRKACB

Image



Western Blot analysis of Jurkat cells using PKA α / β cat Polyclonal Antibody



Western Blot analysis of K562 cells using PKA α/β cat Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.