





# PRKACA/PRKACB Antibody

<b>Product Code</b>	CSB-PA003798
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P17612/P22694
<b>Immunogen</b>	Synthesized peptide derived from the N-terminal region of Human PKA $\alpha$ / $\beta$ cat.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	PRKACA; PKACA; cAMP-dependent protein kinase catalytic subunit alpha; PKA C-alpha; PRKACB; cAMP-dependent protein kinase catalytic subunit beta; PKA C-beta
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	PRKACA/PRKACB

## Image



Western Blot analysis of Jurkat cells using PKA $\alpha$ / $\beta$  cat Polyclonal Antibody



Western Blot analysis of K562 cells using  
PKAα/β cat Polyclonal Antibody