





PRKCD Antibody

Product Code	CSB-PA003803
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05655
Immunogen	Synthesized peptide derived from Human PKC δ around the non-phosphorylation site of T507.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PRKCD; Protein kinase C delta type; Tyrosine-protein kinase PRKCD; nPKC-delta
Immunogen Species	Homo sapiens (Human)
Target Names	PRKCD

Image



Western Blot analysis of 3T3 cells using PKC δ Polyclonal Antibody



Western Blot analysis of various cells using PKC δ Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.