





PRKCZ Antibody

Product Code	CSB-PA003805
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05513
Immunogen	Synthesized peptide derived from Human PKC ζ around the non-phosphorylation site of T560.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PRKCZ; PKC2; Protein kinase C zeta type; nPKC-zeta
Immunogen Species	Homo sapiens (Human)
Target Names	PRKCZ

Image



Western Blot analysis of COS7 cells using PKC ζ Polyclonal Antibody



Western Blot analysis of Hela cells using PKC ζ Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.