





PROKR1 Antibody

Product Code	CSB-PA003812
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCW9
Immunogen	Synthesized peptide derived from the N-terminal region of Human PKR1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PROKR1; GPR73; PKR1; Prokineticin receptor 1; PK-R1; G-protein coupled receptor 73; G-protein coupled receptor ZAQ; GPR73a
Immunogen Species	Homo sapiens (Human)
Target Names	PROKR1

Image



Western Blot analysis of COLO205 cells using PROKR1 Polyclonal Antibody



Western Blot analysis of Jurkat cells using PKR1 Polyclonal Antibody