





PLK5 Antibody

Product Code	CSB-PA003821
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q496M5
Immunogen	Synthesized peptide derived from the Internal region of Human PLK-5.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PLK5; PLK5P; FG060302; Inactive serine/threonine-protein kinase PLK5; Polo-like kinase 5; PLK-5
Immunogen Species	Homo sapiens (Human)
Target Names	PLK5

Image



Western Blot analysis of COLO205 cells using PLK-5 Polyclonal Antibody



Western Blot analysis of Jurkat K562 NIH-3T3 cells using PLK-5 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.