





PMS1 Antibody

Product Code	CSB-PA003824
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54277
Immunogen	Synthesized peptide derived from the Internal region of Human PMS1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	PMS1; PMSL1; PMS1 protein homolog 1; DNA mismatch repair protein PMS1
Immunogen Species	Homo sapiens (Human)
Target Names	PMS1

Image



Western Blot analysis of Jurkat cells using PMS1 Polyclonal Antibody



Western Blot analysis of various cells using PMS1 Polyclonal Antibody

